


Product Name: Recombinant Human UHRF2 Protein
Catalog #: PEH7869



Summary

Name	UHRF2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96PU4
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSTSARARLI DPGFGIYKVN ELVDARDVGL GAWFEAHIHS VTRASDGQSR GKTPLKNGSS CKRTNGNIKH KSKENTNKLD SVPSTSNSDC VAADEDVIYH IQYDEYPESG TLEMNVKDLR PRARTILKWN ELNVGDVVMV NYNVESPGQR GFWFDAEITT LKTISRTKKE LRVKIFLGGS EGTLDNCKII SVDEIFKIER PGAHPLSFAD GKFLRRNDPE CDLCGGDPEK KCHSCSCRVC GGKHEPNMQL LCDECNVAYH IYCLNPPLDK VP EEYWYCP

image



Background

Alternative Names	E3 ubiquitin-protein ligase UHRF2 (EC 2.3.2.27) (Np95/ICBP90-like RING finger protein) (Np95-like RING finger protein) (Nuclear protein 97) (Nuclear zinc finger protein Np97) (RING finger protein 107) (RING-type E3 ubiquitin
--------------------------	--

Product Name: Recombinant Human UHRF2 Protein
Catalog #: PEH7869



transferase UHRF2) (Ubiquitin-like PHD and RING finger domain-containing protein 2) (Ubiquitin-like-containing PHD and RING finger domains protein 2)

Note

For Research Use Only , Not for Diagnostic Use