


Product Name: Recombinant Human GSDMA Protein
Catalog #: PEH7871



Summary

Name	GSDMA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96QA5
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHICNDNMQT FPPGEKSGEE KVILIQASDV GDVHEGFRTL KEEVQRETQQ VEKLSRVGQS SLLSSLSKLL GKKKELQDLE LALEGALDKG HEVTLEALPK DVLLSKEAVG AILYFVGALT ELSEAQQKLL VKSMEKKILP VQLKLVESTM EQNFLLDKEG VFPLQPELLS SLGDEELTLT

image



Background

Alternative Names Gasdermin-A (Gasdermin-1) [Cleaved into: Gasdermin-A, N-terminal (GSDMA-NT); Gasdermin-A, C-terminal (GSDMA-CT)]

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GSDMA Protein
Catalog #: PEH7871

