


Product Name: Recombinant Human ALPK1 Protein
Catalog #: PEH7873



Summary

Name	ALPK1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96QP1
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CSTALSEELE NDREGRAMHS LHSQLHDLSL QEPNNDNLEP SQNQPPQQMP LTPFSPHNTG GFLAPGAGL LEGAPEGIQE VRNMGPRNTS AHSRPSYRSA SWSSDSGRPK NMGTHPSVQK EEAFEIIVEF PETNCDVKDR QGKEQGEEIS ERGAGPTFKA

image



Background

Alternative Names	Alpha-protein kinase 1 (EC 2.7.11.1) (Chromosome 4 kinase) (Lymphocyte alpha-protein kinase)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human ALPK1 Protein
Catalog #: PEH7873

