


Product Name: Recombinant Human PEBP4 Protein
Catalog #: PEH7875



Summary

Name	PEBP4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96S96
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PELGNIGCKV VPDCNNYRQK ITSWMEPIVK FPGAVDGATY ILVMVDPDAP SRAEPRQRFW RHWLVTDIKG ADLKKGKIQG QELSAYQAPS PPAHSGFHRY QFFVYLQEGK VISLLPKENK TRGSWKMDRF LNRFHLEPE ASTQFMTQNY QDSPTLQAPR ERASEPKHKN

image



Background

Alternative Names Phosphatidylethanolamine-binding protein 4 (PEBP-4) (hPEBP4) (Protein cousin-of-RKIP 1)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PEBP4 Protein
Catalog #: PEH7875

