



Product Name: Recombinant Human DCLRE1C Protein
Catalog #: PEH7876



Summary

Name	DCR1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96SD1
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VVEILKPLCR SSQSTEPKYK PLGKLKRART VHRDSEEDD YLFDDPLPIP LRHKVPYPET FHPEVFSMTA VSEKQPEKLR QTPGCCRAEC MQSSRFTNFV DCEESNSESE EEVGIPASLQ GDLGSVLHLQ KADGDVPQWE VFFKRNDIT DESLENFPSS TVAGGSQSPK LFSDDSGEST HISSQNSSQS THITEQGSQG WDSQSDTVLL SSQERNSGDI TSLDKADYRP TIKENIPASL MEQNVICPKD TYSDLKSRDK DVTIVPSTGE PTTLSSETHI PEEKSLLNLS TNADSQSSSD FEVPSTPEAE LPKREHLQYL YEKLATGESI AVKKRKCSLL

image



Background

Alternative Names	Protein artemis (EC 3.1.-.-) (DNA cross-link repair 1C protein) (Protein A-SCID) (SNM1 homolog C) (hSNM1C) (SNM1-like protein)
--------------------------	--

Product Name: Recombinant Human DCLRE1C Protein
Catalog #: PEH7876



Note

For Research Use Only , Not for Diagnostic Use