


Product Name: Recombinant Human PCYT2 Protein
Catalog #: PEH7877



Summary

Name	PCY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99447
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MIRNGRGAAG GAEQPGPGGR RAVRVWCDGC YDMVHYGHSN QLRQARAMGD YLIVGVHTDE EIAKHKGPPV FTQEERYKMV QAIKWVDEVV PAAPYVTTLT TLDKYNCDFC VHGNITLTV DGRDYEYEVK QAGRYRECKR TQGVSTDLV GRMLLVTKAH HSSQEMSSEY REYADSFSGK PGGRNPWTGV SQFLQTSQKI IQFASGKEPQ PGETVIYVAG AFDLFHIGHV DFLEKVHRLA ERPYIIAGLH FDQEVNHYKG KNPIMNLHE RTLSVLACRY VSEVVIGAPY AVTAELLSHF KVDLVCHGKT EIIPDRDGS PYQEPKRRGI FRQIDSGSNL TTDLIVQRRI TNRLEYEARN QKKEAKELAF LEARQQAAQ PLGERDGDF

image



Background

Alternative Names	Ethanolamine-phosphate	cytidyltransferase	(EC	2.7.7.14)
--------------------------	------------------------	--------------------	-----	-----------

Product Name: Recombinant Human PCYT2 Protein
Catalog #: PEH7877



(CTP:phosphoethanolamine cytidyltransferase) (Phosphorylethanolamine transferase)

Note

For Research Use Only , Not for Diagnostic Use