


Product Name: Recombinant Human RIT2 Protein
Catalog #: PEH7881



Summary

Name	RIT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99578
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEVENEASCS PGSASGGSRE YKVVMLGAGG VGKSAMTMQF ISHQFPDYHD PTIEDAYKTQ VRIDNEPAYL DILDTAGQAE FTAMREQYMR GGEGFIICYS VTDRQSFQEA AKFKELIFQV RHTYEIPLVL VGNKIDLEQF RQVSTEEGLS LAQEYNCGFF ETSALRFCI DDAFHGLVRE IRKKESMPSL MEKCLKRKDS LWKKLKGSLK KKRENMT

image



Background

Alternative Names GTP-binding protein Rit2 (EC 3.6.5.2) (Ras-like protein expressed in neurons) (Ras-like without CAAX protein 2)

Note

Product Name: Recombinant Human RIT2 Protein
Catalog #: PEH7881



For Research Use Only , Not for Diagnostic Use