


Product Name: Recombinant Human FEV Protein
Catalog #: PEH7882



Summary

Name	FEV_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99581
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LINMYLPDPV GDGLFKDGKN PSWGPLSPAV QKGSQIQLW QFLELLADR ANAGCIAWEG GHGEFKLTDP DEVARRWGER KSKPNMNYDK LSRALRYYYD KNIMSKVHGK RYAYRFDQFQ LAQACQPPPA HAHAAAAAAAAA AAAAAQDGAL

image



Background

Alternative Names Protein FEV (Fifth Ewing variant protein) (PC12 ETS domain-containing transcription factor 1) (PC12 ETS factor 1) (Pet-1)

Note

For Research Use Only , Not for Diagnostic Use