


Summary

Name	SCAFB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99590
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TSGTSNTRGS RRPAMTTPT RRSTRNTRAE TASQSQRSPI SDNSGCDAPG NSNP SLSVPS SAESEKQTRQ APKRKSVRRG RKPPLLKKKL RSSVAAPEKS SSNDSVDEET AESDTSPVLE KEHQPDVDSS NICTVQTHVE NQSANCLKSC NEQIEESEKH TANYDTEERV GSSSESCAQ DLPVLVGEEG EVKKLENTGI EANVLCLESE ISENILEKGG DPLEKQDQIS GLSQSEVKTD VCTVHLPNDF PTCLTSESKV YQPVSCPLSD LSENVESVVN EEKITESSLV EITEHKDFTL KTEELIESPK

image



Background

Alternative Names	Protein SCAF11 (CTD-associated SR protein 11) (Renal carcinoma antigen NY-REN-40) (SC35-interacting protein 1) (SR-related and CTD-associated factor 11) (SRSF2-interacting protein) (Serine/arginine-rich splicing factor 2-
--------------------------	---

Product Name: Recombinant Human SCAF11 Protein
Catalog #: PEH7883



interacting protein) (Splicing factor, arginine/serine-rich 2-interacting protein) (Splicing regulatory protein 129) (SRp129)

Note

For Research Use Only , Not for Diagnostic Use