


Product Name: Recombinant Human PTPN18 Protein
Catalog #: PEH7887



Summary

Name	PTN18_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & N-His
Accession #	Q99952
Host	E.coli
Species	Human
Molecular Mass	64 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSRSLDSARS FLERLEARGG REGAVLAGEF SDIQACSAAW KADGVCSTVA GSRPENVRKN RYKDVLPYDQ TRVILSLLQE EGHSDYINGN FIRGVDGSLA YIATQGGLPH TLLDFWRLVW EFGVKVILMA CREIENGRKR CERYWAQEQE PLQTGLFCIT LIKEKWLNE IMLRTLKVTF QKESRSVYQL QYMSWPDRGV PSSPDHMLAM VEEARRLQGS GPEPLCVHCS AGCGRTGVLC TVDYVRQLLL TQMIPPDFSL FDVVLKMRKQ RPAAVQTEEQ YRFLYHTVAQ MFCSTLQNAS PHYQNIKENC APLYDDALFL RTPQALLAIP RPPGGVLRSI SVPGSPGHAM ADTYAVVQKR GAPAGAGSGT QTGTGTGTGA RSAEEAPLYS KVTBRAQRP AHAEDARGTL PGRVPADQSP AGSGAYEDVA GGAQTGGLGF NLRIGRPKGP RDPPAEWTRV

image



Background

Product Name: Recombinant Human PTPN18 Protein
Catalog #: PEH7887



Alternative Names

Tyrosine-protein phosphatase non-receptor type 18 (EC 3.1.3.48) (Brain-derived phosphatase)

Note

For Research Use Only , Not for Diagnostic Use