


Product Name: Recombinant Human SH3GL1 Protein
Catalog #: PEH7889



Summary

Name	SH3G1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99961
Host	E.coli
Species	Human
Molecular Mass	42 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSVAGLKKQF YKASQLVSEK VGGAEGTKLD DDFKEMEKKV DVTSKAVTEV LARTIEYLQP NPASRAKLTM LNTVSKIRGQ VKNPGYPQSE GLLGECMIRH GKELGGESNF GDALLDAGES MKRLAEVKDS LDIEVKQNF I DPLQNLCEKD LKEIQHHLKK LEGRRLDYFDY KKKRQGGKIPD EELRQALEKF EESKEVAETS MHNLETDIE QVSQLSALVD AQLDYHRQAV QILDELAEKL KRRMREASSR PKREYKPKPR EPFDLGEPEQ SNGGFPCTTA PKIAASSSFR SSDKPIRTPS RSMPPLDQPS CKALYDFEPE NDGELGFHEG DVITLTNQID ENWYEGMLDG QSGFFPLSYV EVLVPLPQ

image



Background

Alternative Names Endophilin-A2 (EEN fusion partner of MLL) (Endophilin-2) (Extra eleven-

Product Name: Recombinant Human SH3GL1 Protein
Catalog #: PEH7889



nineteen leukemia fusion gene protein) (EEN) (SH3 domain protein 2B) (SH3 domain-containing GRB2-like protein 1)

Note

For Research Use Only , Not for Diagnostic Use