


Product Name: Recombinant Human FERMT1 Protein
Catalog #: PEH7893



Summary

Name	FERM1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BQL6
Host	E.coli
Species	Human
Molecular Mass	66 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LPNLKMOVRLR VSFSAVVFKA VSDICKILNI RRSEELSLK PSGDYFKKKK KKDKNNKEPI IEDILNLESS PTASGSSVSP GLYSKTMTPI YDPINGTPAS STMTWFSDSP LTEQNCSILA FSQPPQSPEA LADMYQPRSL VDKAKLNAGW LDSSRSLMEQ GIQEDEQLLL RFKYYSFFDL NPKYDAVRIN QLYEQARWAI LLEEIDCTEE EMLIFAALQY HISKLSLSAE TQDFAGESEV DEIEAALSNI EVTLEGGKAD SLLEDITDIP KLADNLKLFR PKKLLPKAFK QYWFIFKDT IAYFKNKELE QGEPLEKLNL RGCEVVPDVN VAGRKFGIKL LIPVADGMNE MYLRCDHENQ YAQWMAACML ASKGKTMADS SYQPEVLNLI SFLRMKNRNS ASQVASSLEN MDMNPECFVS PRCAKRHKSK QLAARILEAH QNVAQMPLVE AKLRFIQAWQ SLPEFGLTY LVRFKGSKKD DILGVSYNRL IKIDAATGIP VTTWRFTNIK QWNVNWETRQ VVIEFDQNVF TAFTCLSADC KIVHEYIGGY IFLSTRSKDQ

image



Product Name: Recombinant Human FERMT1 Protein
Catalog #: PEH7893



Background

Alternative Names

Fermitin family homolog 1 (Kindlerin) (Kindlin syndrome protein) (Kindlin-1)
(Unc-112-related protein 1)

Note

For Research Use Only , Not for Diagnostic Use