


Product Name: Recombinant Human NAT9 Protein
Catalog #: PEH7898



Summary

Name	NAT9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BTE0
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MRLNQNTLLL GKKVVLVPYT SEHVPSRYHE WMKSEELQRL TASEPLTLEQ EYAMQCSWQE DADKCTFIVL DAEKWQAQPG ATEESCMVGD VNLFLTDLED LTLGEIEVMI AEPSCRGKGL GTEAVLAML S YGVTTLGLTK FEAIGQGNE PSIRMFQKLH FEQVATSSVF QEVTLRLTVS ESEHQWLLEQ TSHVEEKPYR DGSAEPC

image



Background

Alternative Names N-acetyltransferase 9 (EC 2.3.1.-) (Embryo brain-specific protein)

Note

Product Name: Recombinant Human NAT9 Protein
Catalog #: PEH7898



For Research Use Only , Not for Diagnostic Use