


Product Name: Recombinant Human PAXX Protein
Catalog #: PEH7899



Summary

Name	PAXX_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BUH6
Host	E.coli
Species	Human
Molecular Mass	22 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MDPLSPPLCT LPPGPEPPRF VCYCEGEESG EGDRGGFNLY VTDAAELWST CFTPDSLAL KARFGLSAAE DITPRFRAAC EQQAVALTQ EDRASLTLSG GPSALAFDLS KVPGPAAAPR LRALTGLAK RVWSLERRLA AAETAVSPR KSPRPAGPQL FLDPDPQQRG GPGPGVRRRC PGESLINPGF KSKKPAGGVD FDET

image



Background

Alternative Names Protein PAXX (Paralog of XRCC4 and XLF) (XRCC4-like small protein)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PAXX Protein
Catalog #: PEH7899

