


Product Name: Recombinant Human TRIM62 Protein
Catalog #: PEH7902



Summary

Name	TRI62_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BVG3
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EQHQVTGIDD AFDELQRELK DQLQALQDSE REHTEALQLL KRQLAETKSS TKSLRTTIGE AFERLHRLLR ERQKAMLEEL EADTARTLTD IEQKVQRYSQ QLRKVQEGAQ ILQERLAETD RHTFLAGVAS LSERLKGKIH

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRIM62 (EC 2.3.2.27) (RING-type E3 ubiquitin transferase TRIM62) (Tripartite motif-containing protein 62)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use