


Product Name: Recombinant Human PBDC1 Protein
Catalog #: PEH7903



Summary

Name	PBDC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BVG4
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAATSGTDEP VSGELVSAH ALSLPAESYG NDPDIEMAWA MRAMQHAENV YKLISSVDPQ FLKLTQVDDQ IYSEFRKNFE TLRIDVLDPE ELKSESAKEK WRPFCLKFNG IVEDFNYGTL LRLDCSQGYT EENTIFAPRI QFFAIEIARN REGYNKAVYI SVQDKGEKEG VNNGGEKRAD SGEEENTKNG GEKGADSGEE KEEGINREDK TDKGGEKGKE ADKEINKSGE KAM

image



Background

Alternative Names	Protein PBDC1 (Polysaccharide biosynthesis domain-containing protein 1)
--------------------------	---

Note

Product Name: Recombinant Human PBDC1 Protein
Catalog #: PEH7903



For Research Use Only , Not for Diagnostic Use