


Product Name: Recombinant Human TRAIP Protein
Catalog #: PEH7904



Summary

Name	TRAIP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BWF2
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ENVLDAEFLK NELDNVRAQL SQKDKEKRDS QVIIDTLRDT LEERNATVVS LQQALGKAEM LCSTLKKQMK YLEQQQDETQ QAQEEARRLR SKMKTMEQIE LLLQSRPEV EEMIRDMGVG QSAVEQLAVY CVSLKKEYEN LKEARKASGE VADKLRKDLF SSRSLQTVY SELDQAKLEL KSAQKDLQSA

image



Background

Alternative Names	E3 ubiquitin-protein ligase TRAIP (EC 2.3.2.27) (RING finger protein 206) (TRAF-interacting protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TRAIP Protein
Catalog #: PEH7904

