


Product Name: Recombinant Human SLC26A6 Protein
Catalog #: PEH7911



Summary

Name	S26A6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BXS9
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PHYSVLGQVP DTDIYRDVAE YSEAKEVRGV KVFRSSATVY FANAEFYSDA LKQRCGVDVD FLISQKKLL KKQEQLKLKQ LQKEEKLRLKQ AASPKGASVS INVNTSLEDM RSNNVEDCKM MQVSSGDKME DATANGQEDS KAPDGSTLKA LGLPQPDFHS LILDLGALSF VDTVCLKSLK NIFHDFREIE VEVYMAACHS PVVSQLEAGH FFDASITKKH LFASVHDAVT FALQHPRPVP

image



Background

Alternative Names	Solute carrier family 26 member 6 (Anion exchange transporter) (Pendrin-like protein 1) (Pendrin-L1)
--------------------------	--

Note

Product Name: Recombinant Human SLC26A6 Protein
Catalog #: PEH7911



For Research Use Only , Not for Diagnostic Use