


**Product Name: Recombinant Human PPP1R12C Protein**  
**Catalog #: PEH7915**



## Summary

|                              |                                                                                                                                                                                                                |
|------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Name</b>                  | PP12C_HUMAN                                                                                                                                                                                                    |
| <b>Purity</b>                | Greater than 85% as determined by reducing SDS-PAGE                                                                                                                                                            |
| <b>Endotoxin level</b>       | Please contact with the lab for this information                                                                                                                                                               |
| <b>Tag</b>                   | C-His                                                                                                                                                                                                          |
| <b>Accession #</b>           | Q9BZL4                                                                                                                                                                                                         |
| <b>Host</b>                  | E.coli                                                                                                                                                                                                         |
| <b>Species</b>               | Human                                                                                                                                                                                                          |
| <b>Molecular Mass</b>        | 18 kDa                                                                                                                                                                                                         |
| <b>Buffer</b>                | Phosphate buffered saline                                                                                                                                                                                      |
| <b>Shipping</b>              | The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.                                                                                             |
| <b>Stability&amp;Storage</b> | Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles. |
| <b>Sequence</b>              | EWRGPAEGEE AEPADRSQES STLEGGPSAR RQRWQRDLNP EPEPESEEPD<br>GGFRTLAEEL RRENERLREA LTETTLRLAQ LKVELERATQ RQERFAERPA<br>LLELERFERR ALERKAAELE EELKALSDLR ADNQRKLDEN                                                |

## image



## Background

|                          |                                                                                                                                                          |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Alternative Names</b> | Protein phosphatase 1 regulatory subunit 12C (Protein phosphatase 1 myosin-binding subunit of 85 kDa) (Protein phosphatase 1 myosin-binding subunit p85) |
|--------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|

## Note

For Research Use Only , Not for Diagnostic Use

**Product Name: Recombinant Human PPP1R12C Protein**  
**Catalog #: PEH7915**

---

