


Product Name: Recombinant Human YIPF3 Protein
Catalog #: PEH7920



Summary

Name	YIPF3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZM5
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFEENIQGGG SAVIDMENMD DTSGSSFEDM GELHQRLREE EVDADAADAA AAEEEDGEFL GMKGFKGQLS RQVADQMWQA

image



Background

Alternative Names	Protein YIPF3 (Killer lineage protein 1) (Natural killer cell-specific antigen KLIP1) (YIP1 family member 3) [Cleaved into: Protein YIPF3, 36 kDa form III]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use