


Product Name: Recombinant Human C19orf33 Protein
Catalog #: PEH7922



Summary

Name	IMUP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZP8
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MEFDLGAAL PTSQKPGVGA GHGGDPKLSP HKVQGRSEAG AGPGPKQGHH SSSDSSSSSS DSDTDVKSHA AGSKQHESIP GKAKKPKVKK KEKGKKEKGK KKEAPH

image



Background

Alternative Names	Immortalization up-regulated protein (Hepatocyte growth factor activator inhibitor type 2-related small protein) (H2RSP) (HAI-2-related small protein)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use