


Product Name: Recombinant Human NIF3L1 Protein
Catalog #: PEH7925



Summary

Name	NIF3L_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9GZT8
Host	E.coli
Species	Human
Molecular Mass	58 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLVIRALENR VGIYSPHTAY DAAPQGVNNW LAKGLGACTS RPIHPSKAPN YPTEGNHRVE FNVNYTQDLD KVMSAVKGID GVSVTSFSAR TGNEEQTRIN LNCTQKALMQ VVDFLSRNKQ LYQKTEILSL

image



Background

Alternative Names	NIF3-like protein 1 (Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 1 protein)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use