


Summary

Name	PAIP1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H074
Host	E.coli
Species	Human
Molecular Mass	54 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSDGFDRAPG AGRGRSRGLG RGGGGPEGGG FPNGAGPAER ARHQPPQPKA PGFLQPPPLR QPRTTPPPGA QCEVPASPQR PSRPGALPEQ TRPLRAPPSS QDKIPQQNSE SAMAKPQVVV APVLMKLSV NAPEFYPSGY SSSYTESYED GCEDYPTLSE YVQDFLNHLT EQPGSFETEI EQFAETLNGC VTTDDALQEL VELIYQQATS IPNFSYMGAR LCNYLSHHLT ISPQSGNFRQ LLLQRCRTEY EVKDQAAKGD EVTRKRFHAF VLFLGELYLN LEIKGTNGQV TRADILQVGL RELLNALFSN PMDDNLICAV KLLKLTGSVL EDAWKEKGKM DMEEIIQRIE NVVLDANCSR DVKQMLLKLK ELRSSNWGRV HATSTYREAT PENDPNYFMN EPTFYTSDGV PFTAADPDYQ EKYQELLERE DFFPDYEENG TDLSGAGDPY LDDIDDEMDP EIEEAYEKFC LESERKRKQ

image



Background

Product Name: Recombinant Human PAIP1 Protein
Catalog #: PEH7927



Alternative Names

Polyadenylate-binding protein-interacting protein 1 (PABP-interacting protein 1) (PAIP-1) (Poly(A)-binding protein-interacting protein 1)

Note

For Research Use Only , Not for Diagnostic Use