


Product Name: Recombinant Human PCBD2 Protein
Catalog #: PEH7928



Summary

Name	PHS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H0N5
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NQAILDLKAA GWSELSERDA IYKEFSFHNF NQAFGFMSRV ALQAEKMNH PEWFNVYNKV QITLTSHDCG ELTKKDVKLA

image



Background

Alternative Names	Pterin-4-alpha-carbinolamine dehydratase 2 (PHS 2) (EC 4.2.1.96) (4-alpha-hydroxy-tetrahydropterin dehydratase 2) (DcoH-like protein DCoHm) (Dimerization cofactor of hepatocyte nuclear factor 1 from muscle) (HNF-1-alpha dimerization cofactor)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human PCBD2 Protein
Catalog #: PEH7928

