


Product Name: Recombinant Human ATG5 Protein
Catalog #: PEH7929



Summary

Name	ATG5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1Y0
Host	E.coli
Species	Human
Molecular Mass	33 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTDDKDVLRD VWFGRIPTCF TLYQDEITER EAEPYLLLP RVSYLTLVTD KVKKHFQKVM RQEDISEIWF EYEGTPLKWH YPIGLLFDLL ASSSALPWN TVHFKSFPEK DLLHCPSKDA IEAHFMSCMK EADALKHKSQ VINEMQKKDH KQLWMGLQND RFDQFWAINR KLMEYPAEEN GFRYIPFRIY QTTTERPFIQ KLFRPVAADG QLHTLGDLLK EVCPSAIDPE DGEKKNQVMI HGIEPMLETP LQWLSEHLSY PDNFLHISII PQPTD

image



Background

Alternative Names Autophagy protein 5 (APG5-like) (Apoptosis-specific protein)

Product Name: Recombinant Human ATG5 Protein
Catalog #: PEH7929



Note

For Research Use Only , Not for Diagnostic Use