


Product Name: Recombinant Human SH3BGRL3 Protein
Catalog #: PEH7931



Summary

Name	SH3L3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H299
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSGRLRVYSTS VTGSREIKSQ QSEVTRILDG KRIQYQLVDI SQDNALRDEM RALAGNPKAT PPQIVNGDQY CGDYELFVEA VEQNTLQEFL KLA

image



Background

Alternative Names SH3 domain-binding glutamic acid-rich-like protein 3 (SH3 domain-binding protein 1) (SH3BP-1)

Note

For Research Use Only , Not for Diagnostic Use