



Product Name: Recombinant Human SLK Protein
Catalog #: PEH7932



Summary

Name	SLK_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2G2
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RRFEQEMMSK KRQYDQEIEN LEKQQKQTIE RLEQEHTNRL RDEAKRIKGE QEKELSKFQN MLKNRKKEVI NEVEKAPKEL RKELMKRRKE ELAQSQHAQE QEFVQKQQQE LDGSLKKIIQ QQKAELANIE RECLNKKQQL MRAREAAIWE LEERHLQEKH QLLKQQLKDQ YFMQRHQLLK RHEKETEQQM RYNQRLIEEL KNRQTQERAR LPKIQRSEAK TRMAMFKKSL RINSTATPDQ DRDKIKQFAA QEEKRQKNER MAQHQKHENQ MRDLQLQCEA NVRELHQLQN EKCHLLVEHE

image



Background

Alternative Names	STE20-like serine/threonine-protein kinase (STE20-like kinase) (hSLK) (EC 2.7.11.1) (CTCL tumor antigen se20-9) (STE20-related serine/threonine-protein kinase) (STE20-related kinase) (Serine/threonine-protein kinase 2)
--------------------------	--

Product Name: Recombinant Human SLK Protein
Catalog #: PEH7932



Note

For Research Use Only , Not for Diagnostic Use