


Product Name: Recombinant Human RAB3GAP2 Protein
Catalog #: PEH7934



Summary

Name	RBGPR_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H2M9
Host	E.coli
Species	Human
Molecular Mass	12 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSQELESVDE GLLQFCANKL KLLQLYESVS QLNSLDFHLD TPFSDNDLAL LLRLDEKELL KLQALLEKYK QENTRTNVRF

image



Background

Alternative Names	Rab3 GTPase-activating protein non-catalytic subunit (RGAP-iso) (Rab3 GTPase-activating protein 150 kDa subunit) (Rab3-GAP p150) (Rab3-GAP150) (Rab3-GAP regulatory subunit)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use