


Product Name: Recombinant Human CHMP4B Protein
Catalog #: PEH7935



Summary

Name	CHM4B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H444
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MSVFGKLFGA GGGKAGKGGP TPQEAIQRLR DTEEMLSKKQ EFLEKKIEQE LTAACKHGTK NKRAALQALK RKKRYEKQLA QIDGTLSTIE FQREALENAN TNTEVLKNMG YAAKAMKAAH DNMDIDKVDE LMQDIADQQE LAEEISTAIS KPVGFGEEDF EDELMAELEE LEQEELDKNL LEISGPETVP LPNVPSIALP SKPAKKKEEE DDDMKELENW AGSM

image



Background

Alternative Names	Charged multivesicular body protein 4b (Chromatin-modifying protein 4b) (CHMP4b) (SNF7 homolog associated with Alix 1) (SNF7-2) (hSnf7-2) (Vacuolar protein sorting-associated protein 32-2) (Vps32-2) (hVps32-2)
--------------------------	---

Product Name: Recombinant Human CHMP4B Protein
Catalog #: PEH7935



Note

For Research Use Only , Not for Diagnostic Use