


Product Name: Recombinant Human SENP3 Protein
Catalog #: PEH7937



Summary

Name	SENP3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H4L4
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SNVCSIGDHV AQELFQGS DL GMAEEAERPG EKAGQHSPLR EEHVTCVQSI LDEFLQTYGS LIPLSTDEVV EKLEDIFQQE

image



Background

Alternative Names Sentrin-specific protease 3 (EC 3.4.22.-) (SUMO-1-specific protease 3)
(Sentrin/SUMO-specific protease SENP3)

Note

For Research Use Only , Not for Diagnostic Use