


## Summary

<b>Name</b>	PPIP2_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9H939
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	40 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	MTRSLFKGNF WSADILSTIG YDNIIQHLNN GRKNCKEFED FLKERA AIEE RYGKDLLNLS RKKPCGQSEI NTLKRALEVF KQQVDNVAQC HIQLAQLSRE EARKMEEFRE KQKLQRKKTE LIMDAIHKQK SLQFKKTMDA KKNYEQKCRD KDEAEQAVSR SANLVNPKQQ EKLFVKLATS KTAVEDSDKA YMLHIGTLDK VREEWQSEHI KACEAFEAEQ CERINFFRNA LWLHVNQLSQ QCVTSDEMYE QVRKSLEMCS IQRDIEYFVN QRKTGQIPPA PIMYENFYSS QKNAVPAGKA TGPNLARRGP LPIPKSSPDD PNYSLVDDYS LLYQ

## image



## Background

<b>Alternative Names</b>	Proline-serine-threonine phosphatase-interacting protein 2 (PEST phosphatase-interacting protein 2)
--------------------------	---

**Product Name: Recombinant Human PSTPIP2 Protein**  
**Catalog #: PEH7939**

---



**Note**

For Research Use Only , Not for Diagnostic Use