


Product Name: Recombinant Human AGO3 Protein
Catalog #: PEH7942



Summary

Name	AGO3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H9G7
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KLLANCFQVE IPKIDVLYE VDIKPDKCPR RVNREVVDSM VQHFKVTIFG DRRPVYDGKR SLYTANPLPV ATTGVDLDTV LPGEGGKDRP FKVSIKFSR VSWHLLHEVL TGRTLPEPLE

image



Background

Alternative Names	Protein argonaute-3 (Argonaute3) (hAgo3) (EC 3.1.26.n2) (Argonaute RISC catalytic component 3) (Eukaryotic translation initiation factor 2C 3) (eIF-2C 3) (eIF2C 3)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human AGO3 Protein
Catalog #: PEH7942

