


Product Name: Recombinant Human PPCS Protein
Catalog #: PEH7944



Summary

Name	PPCS_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HAB8
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAEMDPVAEF PQPPGAARWA EVMARFAARL GAQGRRVVLV TSGGTKVPLE ARPVRFLDNF SSGRRGATSA EAFLAAGYGV LFLYRARSFAF PYAHRFPQQT WLSALRPSGP ALSGLLSLEA EENALPGFAE ALRSYQEAAA AGTFLAVEFT TLADYLHLLQ AAAQALNPLG PSAMFYLAAA VSDFYVPVSE MPEHKIQSSG GPLQITMKMV PKLLSPLVKD WAPKAFIISF KLETDPAIVI NRARKALEIY QHQQVVANIL ESRQSFVIV TKDSETKLLL SEEEIEKGVE IEEKIVDNLQ SRHTAFIGDR N

image



Background

Alternative Names Phosphopantothenate--cysteine ligase (EC 6.3.2.51)
(Phosphopantothenoylcysteine synthetase) (PPC synthetase)

Product Name: Recombinant Human PPCS Protein
Catalog #: PEH7944



Note

For Research Use Only , Not for Diagnostic Use