

Product Name: Recombinant Human ZBTB20 Protein
Catalog #: PEH7949



Summary

Name	ZBT20_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HC78
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TPRGTPESGT SGQSSDTESG YLQSH PQHSV DRIYSALYAC SMQNGSGERS FYSGAVVSHH ETALGLPRDH H MEDPSWITR IHERSQQMER YLSTTPETTH CRKQPRPVRI QTLVGNIIHK QEMEDDYDYY GQQRVQILER NESEECTEDT DQAEGTESEP KGESFDSGVS SSIGTEPDSV EQQFGPGAAR DSQAEPQPE QAAEAPAEGG PQTNQLETGA SSPERSNEVE MDSTVITVSN SSDKSVLQQP

image



Background

Alternative Names Zinc finger and BTB domain-containing protein 20 (Dendritic-derived BTB/POZ zinc finger protein) (Zinc finger protein 288)

Note

Product Name: Recombinant Human ZBTB20 Protein
Catalog #: PEH7949



For Research Use Only , Not for Diagnostic Use