


Product Name: Recombinant Human NCOA5 Protein
Catalog #: PEH7950



Summary

Name	NCOA5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCD5
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	CREKEREIIA RQAAKMADEA ILQERERGGP EEGVRGGHPP AIQSLINLLA DNRYLTAET DKIIINYLRER KERLMRSSTD SLPGPISRQP LGATSGASLK TQPSSQPLQS GQVLPSATPT PSAPPTSQQE LQAKILSLFN SGTVTANSSS ASPSVAAGNT PNQNFSTAAN SQPQRSQAS GNQPPSILGQ GGSAQNMGPR

image



Background

Alternative Names Nuclear receptor coactivator 5 (NCoA-5) (Coactivator independent of AF-2) (CIA)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human NCOA5 Protein
Catalog #: PEH7950

