


Product Name: Recombinant Human SMURF1 Protein
Catalog #: PEH7951



Summary

Name	SMUF1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9HCE7
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSIKIRLTVL CAKNLAKKDF FRLPDPFAKI VVDGSGQCHS TDTVKNTLDP KWNQHYDLYV GKTDSITISV WNHKKIHKKQ GAGFLGCVRL LSNAISRLKD TGYQRDLCK LNPSDTDAVR GQIVVSLQTR DRIGTGGSVV DCRGLLENEG TVYEDSGPGR PLSCFMEEPA PYTDSTGAAA GGGNCRFVES PSQDQRLQAQ

image



Background

Alternative Names	E3 ubiquitin-protein ligase SMURF1 (hSMURF1) (EC 2.3.2.26) (HECT-type E3 ubiquitin transferase SMURF1) (SMAD ubiquitination regulatory factor 1) (SMAD-specific E3 ubiquitin-protein ligase 1)
--------------------------	--

Note

Product Name: Recombinant Human SMURF1 Protein
Catalog #: PEH7951



For Research Use Only , Not for Diagnostic Use