


Product Name: Recombinant Human NIT2 Protein
Catalog #: PEH7952



Summary

Name	NIT2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NQR4
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MTSFRLALIQ LQISSIKSDN VTRACSFIRE AATQGAKIVS LPECFNSPYG AKYFPEYAEK IPGESTQKLS EVAKECSIYL IGGSIPEDA GKLYNTCAVF GPDGTLLAKY RKIHLFDIDV PGKITFQESK TLSPGDSFST FDTPYCRVGL GICYDMRFAE LAQIYAQRGC QLLVYPGAFN LTTGPAHWEL LQRSRAVDNQ VYVATASPAR DDKASYVAWG HSTVVNPWGE VLAKAGTEEA IVYSDIDLKK LAEIRQQIPV FRQKRSDLYA VEMKKP

image



Background

Alternative Names Omega-amidase NIT2 (EC 3.5.1.3) (Nitrilase homolog 2)

Product Name: Recombinant Human NIT2 Protein
Catalog #: PEH7952



Note

For Research Use Only , Not for Diagnostic Use