


Product Name: Recombinant Human ACS2 Protein
Catalog #: PEH7953



Summary

Name	ACSA_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NR19
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GPLANGATSV LFEGIPTYPD VNRLWSIVDK YKVTKFYTAP TAIRLLMKFG DEPVTKHSRA SLQVLGTVGE PINPEAWLWY HRVVGQAQRCP IVDTFWQTET GGHMLTPLPG ATPMKPGSAT FPFPGVAPAI LNESGEELEG EAEGYLVFKQ PWPGIMRTVY GNHERFETTY FKKFPGYYVT GDGCQRDQDG YYWITGRIDD MLNVSGHLLS TAEVESALVE HEAVAEAAVV

image



Background

Alternative Names	Acetyl-coenzyme A synthetase, cytoplasmic (EC 6.2.1.1) (Acetate--CoA ligase) (Acetyl-CoA synthetase) (ACS) (AceCS) (Acetyl-CoA synthetase 1) (AceCS1) (Acyl-CoA synthetase short-chain family member 2) (Acyl-activating enzyme) (Propionate--CoA ligase) (EC 6.2.1.17)
--------------------------	---

Product Name: Recombinant Human ACSS2 Protein
Catalog #: PEH7953



Note

For Research Use Only , Not for Diagnostic Use