


Product Name: Recombinant Human POLE3 Protein
Catalog #: PEH7955



Summary

Name	DPOE3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRF9
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MAERPDLNL PNAVITRIK EALPDGVNIS KEARSAISRA ASVFVLYATS CANNFAMKGGK RKTNLASDVL SAMEEMEFQR FVTPLKEALE AYRREQKGKK EASEQKKKDK DKKTDSEEQD KSRDEDNDED EERLEEEEQN EEEVD

image



Background

Alternative Names	DNA polymerase epsilon subunit 3 (Arsenic-transactivated protein) (AsTP) (Chromatin accessibility complex 17 kDa protein) (CHRAC-17) (HuCHRAC17) (DNA polymerase II subunit 3) (DNA polymerase epsilon subunit p17)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human POLE3 Protein
Catalog #: PEH7955

