


Product Name: Recombinant Human DUSP22 Protein
Catalog #: PEH7956



Summary

Name	DUS22_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NRW4
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GLYIGNFKDA RDAEQLSKNK VTHILSVHDS ARPMLEGVKY LCIPAADSPS QNLTRHFKES IKFIHECRLR GESCLVHCLA GVSRSVTLVI AYIMTVTDFG WEDALHTVRA GRSCANPNVG FQRQLQEF EK HEVHQYRQWL KEEYGESPLQ

image



Background

Alternative Names	Dual specificity protein phosphatase 22 (EC 3.1.3.16) (EC 3.1.3.48) (JNK-stimulatory phosphatase-1) (JSP-1) (Low molecular weight dual specificity phosphatase 2) (LMW-DSP2) (Mitogen-activated protein kinase phosphatase x) (MAP kinase phosphatase x) (MKP-x)
--------------------------	--

Note

Product Name: Recombinant Human DUSP22 Protein
Catalog #: PEH7956



For Research Use Only , Not for Diagnostic Use