


Product Name: Recombinant Human LANCL2 Protein
Catalog #: PEH7958



Summary

Name	LANC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-MBP & C-His
Accession #	Q9NS86
Host	E.coli
Species	Human
Molecular Mass	84 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RPPATTDEPG LPFHQDGKII HNFIRRIQTK IKDLLQQMEE GLKTADPHDC SAYTGWTGIA LLYLQLYRVT CDQTYLLRSL DYVKRTLRLN NGRRVTFLCG DAGPLAVGAV IYHKLRSDCE SQECVTKLLQ LQRSVVCQES DLPDELLYGR AGYLYALLYL NTEIGPGTVC ESAIKEVVNA IIESGKTLSE EERKTERCPL LYQWHRKQYV GAAHGMAGIY YMLMQPAAKV DQETLTEMVK PSIDYVRHKK FRSGNYPSSL SNETDRLVHW CHGAPGVIHM LMQAYKVFKE EKYLKEAMEC SDVIWQRGLL RKGYGICHGT AGNGYSFLSL YRLTQDKKYL YRACKFAEWC LDYGAHGCR

image



Background

Alternative Names	LanC-like protein 2 (Testis-specific adriamycin sensitivity protein)
--------------------------	--

Product Name: Recombinant Human LANCL2 Protein
Catalog #: PEH7958



Note

For Research Use Only , Not for Diagnostic Use