


Product Name: Recombinant Human LANCL2 protein
Catalog #: PEH7959



Summary

Name	LANC2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NS86
Host	E.coli
Species	Human
Molecular Mass	27 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHLGGEAEME ERAFVNPFDP YEAAAGALLA SGAAEETGCV RPPATTDEPG LPFHQDGKII HNFIRRIQTK IKDLLQQMEE GLKTADPHDC SAYTGWTGIA LLYLQLYRVT CDQTYLLRSL DYVKRTLRLN NGRRVTFCLG DAGPLAVGAV IYHKLRSDCE SQECVTKLLQ LQRSVVCQES DLPDELLYGR AGYLYALLYL NTEIGPGTVC ESAIKEVVNA

image



Background

Alternative Names LanC-like protein 2 (Testis-specific adriamycin sensitivity protein)

Note

Product Name: Recombinant Human LANCL2 protein
Catalog #: PEH7959



For Research Use Only , Not for Diagnostic Use