


Product Name: Recombinant Human CZIB Protein
Catalog #: PEH7967



Summary

Name	CZIB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NWW4
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MGKIALQLKA TLENITNLRP VGEDFRWYLK MKCGNCGEIS DKWQYIRLMD SVALKGGRGS ASMVQKCKLC ARENSIEILS STIKPYNAED NENFKTIVEF ECRGLEPVDF QPQAGFAAEG VESGTAFSDI NLQEKDWDY DEKAQESVGI YEVTHQFVKC

image



Background

Alternative Names CXXC motif containing zinc binding protein (UPF0587 protein C1orf123)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human CZIB Protein
Catalog #: PEH7967

