


Product Name: Recombinant Human UGGT1 Protein
Catalog #: PEH7971



Summary

Name	UGGG1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NYU2
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GFWDSFKWGF TGQKTEEVKQ DKDDIINIFS VASGHLYERF LRIMMLSVLK NTKTPVKFWF LKNYLSPTFK EFIPYMANEY NFQYELVQYK WPRWLHQQTE KQRIIWGYKI LFLDVLPLV VDKFLFVDAD QIVRTDLKEL RDFNLDGAPY GYTPFCDSRR EMDGYRFWKS GYWASHLAGR KYHISALYVV DLKKFRKIAA GDRLRGQYQG LSQDPNSLSN LDQDLPNMI HQVPIKSLPQ EWLWCETWCD DASKKRAKTI DLCNNPMTKE PKLEAAVRIV PEWQDYDQEI KQLQIRFQKE KETGALYKEK TKEPSREGPQ

image



Background

Alternative Names UDP-glucose:glycoprotein glucosyltransferase 1 (UGT1) (hUGT1) (EC 2.4.1.-)
(UDP--Glc:glycoprotein glucosyltransferase) (UDP-glucose ceramide

Product Name: Recombinant Human UGGT1 Protein
Catalog #: PEH7971



glucosyltransferase-like 1)

Note

For Research Use Only , Not for Diagnostic Use