


Product Name: Recombinant Human MAT2B Protein
Catalog #: PEH7972



Summary

Name	MAT2B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZL9
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VLVTGATGLL GRAVHKEFQQ NNWHAVGCGF RRARPKFEQV NLLDSNAVHH IIHDFQPHVI VHCAERRPD VVENQPDAAS QLNVDASGNL AKEAAAVGAF LIYSSDYVF DGTNPPYREE DIPAPLNLYG KTKLDGEKAV LENNLGAAVL RIPILYGEVE KLEESAVTVM FDKVQFSNKS ANMDHWQQRF PTHVKDVATV CRQLAEKRML DPSIKGTFHW SGNEQMTKYE MACAIADAFN LPSSHLRPIT DSPVLGAQRP RNAQLDCSKL ETLGIGQRTP FRIGIKESLW PFLIDKRWRQ

image



Background

Alternative Names	Methionine adenosyltransferase 2 subunit beta (Methionine adenosyltransferase II beta) (MAT II beta) (Putative dTDP-4-keto-6-deoxy-D-glucose 4-reductase)
--------------------------	---

Product Name: Recombinant Human MAT2B Protein
Catalog #: PEH7972



Note

For Research Use Only , Not for Diagnostic Use