



Product Name: Recombinant Human EHD3 Protein
Catalog #: PEH7973



Summary

Name	EHD3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9NZN3
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LRKLNDLIKR ARLAKVHAYI ISSLKKEMPS VFGKDNKKKE LVNNLAEIYG RIEREHQISP GDFPNLKRMQ DQLQAQDFSK FQPLKSKLLE VVDDMLAHD AQLMVLVRQE ESQRPIQMVK GGAFEGTLHG PFGHGYGEGA GEGIDDAEWV VARDKPMYDE IFYTLSPVDG KITGANAKKE MVRSKLPNSV LGKIWKLADI DKDGMLDDDE FALANHLIKV KLEGHELPNE LPAHLLPPSK

image



Background

Alternative Names	EH domain-containing protein 3 (PAST homolog 3)
--------------------------	---

Note

Product Name: Recombinant Human EHD3 Protein
Catalog #: PEH7973



For Research Use Only , Not for Diagnostic Use