


Product Name: Recombinant Human SLAIN2 Protein
Catalog #: PEH7978



Summary

Name	SLAI2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P270
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DRNPPLSPQS SIDSELSASE LDEDSIGSNY KLNDVTDVQI LARMQEESLR QEYAATTSSR SSGSSCNSTR RGTFSQELD AQLDDEDDN MHHAVYPAVN RFSPSPRNSP RPSPKQSPRN SPRSRSPARG IEYSRVSPQP MISRLQQPRL SLQGHPTDLQ TSNVKNEEKL RRSLPNLSRT SNTQVDSVKS SRSDSNFQVP NGGIPRMQPQ ASAIPSPGKF RSPAAPSPLA LRQPVKAFSN HGSGSPGSQE ITQLTQTSS PGPPMVQSTV SANPPSNINS ATLTRPAGTT AMRSGLPVRPS APSAGGIPVP RSKLAQPVRR

image



Background

Alternative Names	SLAIN motif-containing protein 2
--------------------------	----------------------------------

Product Name: Recombinant Human SLAIN2 Protein
Catalog #: PEH7978



Note

For Research Use Only , Not for Diagnostic Use