


Product Name: Recombinant Human CGN Protein
Catalog #: PEH7981



Summary

Name	CING_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9P2M7
Host	E.coli
Species	Human
Molecular Mass	43 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	RLEEAQRGLA RLGQEQQTLN RALEEEGKQR EVLRRGKAEL EEQKRLLDRT VDRLNKELEK IGEDSKQALQ QLQAQLEDYK EKARREVADA QRQAKDWASE AEKTSGGLSR LQDEIQRRLRQ ALQASQAERD TARLDKELLA QRLQGLEQEA ENKKRSQDDR ARQLKGLEEK VSRLETDELDE EKNTVELLTD RVNRGRDQVD QLRTELMQER SARQDLECDK ISLERQNKDL KTRLASSEGF QKPSASLSQL ESQNQLLQER LQAEEREKTV LQSTNRKLER KVKELSIQIE DERQHVNDQK DQLSLRVKAL KRQVDEAEIE IERLDGLRKK AQRVEVEEQHE VNEQLQARIK

image



Background

Alternative Names	Cingulin
--------------------------	----------

Product Name: Recombinant Human CGN Protein
Catalog #: PEH7981



Note

For Research Use Only , Not for Diagnostic Use