



Summary

Name	EP15R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBC2
Host	E.coli
Species	Human
Molecular Mass	49 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YSLEQDIREK EEAIRQKTSE VQELQNLDLR ETSSLQELEA QKQDAQDRLD EMDQQKAKLR DMLSDVRQKC QDETQMISL KTQIQSQESD LKSQEDDLNR AKSELNRLQQ EETQLEQSIQ AGRVQLETII KSLKSTQDEI NQARSKLSQL HESRQEAHRS LEQYDQVLDG AHGASLTDLA NLSEGVSLAE RGSFGAMDDP FKNKALLFSN NTQELHPDPF QTEDPFKSDP FKGADPFKGD PFQNDPFAEQ QTTSTDPFGG DPFKESDPFR GSATDDFFKK QTKNDPFTSD PFTKNPSLPS KLDPFESSDP FSSSSVSSKG SDPFGTLDPF GSGSFNSAEG FADFSQMSKP PPSGPFTSSL GGAGFSDDPF KSKQDTPALP PKKPAPPRPK PPSGKSTPVS QLGSADFPEA PDPFQPLGAD

image



Background

Product Name: Recombinant Human EPS15L1 Protein
Catalog #: PEH7982



Alternative Names

Epidermal growth factor receptor substrate 15-like 1 (Eps15-related protein)
(Eps15R)

Note

For Research Use Only , Not for Diagnostic Use