


Product Name: Recombinant Human PPARGC1A Protein
Catalog #: PEH7983



Summary

Name	PRGC1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UBK2
Host	E.coli
Species	Human
Molecular Mass	44 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSLLKKLLLA PANTQLSYNE CSGLSTQNHA NHNHRI RTNP AIVKTENSW NKAKSICQQQ KPQRRPCSEL LKYLTTNDDP PHTKPTENRN SSRDKCTSKK KSHTQSQSQH LQAKPTTSL PLTPESPNDP KGSPFENKTI ERTLSVELSG TAGLTPPTTP PHKANQDNPF RASPCLKSSC KTVVPPPSKK PRYSESSGTQ GNNSTKKGPE QSELYAQLSK SSVLTGGHEE RKTKRPSLRL FGDHDYCQSI NSKTEILINI SQELQDSRQL ENKDVSSDWQ GQICSSTDSD QCYLRETLA SKQVSPCSTR KQLQDQEIRA ELNKHFGHPS QAVFDDEADK TGELRDSDFS NEQFSKLPMF INSGLAMDGL FDDSEDESDK

image



Background

Alternative Names	Peroxisome proliferator-activated receptor gamma coactivator 1-alpha
--------------------------	--



Product Name: Recombinant Human PPARGC1A Protein
Catalog #: PEH7983

(PGC-1-alpha) (PPAR-gamma coactivator 1-alpha) (PPARGC-1-alpha) (Ligand effect modulator 6)

Note

For Research Use Only , Not for Diagnostic Use